



Email : sales@laboratoryinstrumentindia.com

Product Name :
Tube Cutter

Product Code :
LBNY-0004-920001



Description :

Tube Cutter, Tube Cutter Equipments, Tube Cutter Tools, Tube Cutter Tool Kits, Tube Cutter Manufacturers, Tube Cutter Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Gabon, Lesotho, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia

Technical Specification :

Tube cutter

There is a quick release mechanism. Used for cutting PVC, copper, brass, aluminum, stainless steel pipes. Adjustable cutting depth design for fast and safe cutting. Good results and precise cutting.

Ergonomic and easy to use, non-slip handle make it comfortable, stable to handle and safe in cutting.

The thickness and accuracy can be adjusted, and the cutting capacity of the tubular cutter is 6-42mm.

The pipe cutter is made of high-strength alloy steel, and the knob contains a spare cutting wheel.

Heavy duty tube cutter

Aluminum material

Suitable for cutting: Copper, Aluminum, Thin stainless steel and PVC

Cutting size :6-42mm OD



Laboratory Instrument India