



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pipe Tube Cutter

Product Code :
LBNY-0004-920002



Description :

Pipe Tube Cutter

Technical Specification :

Tube cutter

Excellent material, strong and durable.

Small design, easy to carry.

Fine workmanship, no burr.

It is suitable for operation in narrow space, cutting copper and aluminum pipes

Mini tube cutter

Aluminum material

Suitable for cutting: Copper

Cutting size: 3-16mm OD



Laboratory Instrument India