



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Hand Drill

**Product Code :**  
LBNY-00011-1010003



**Description :**

It is used for hatching eggs of chicken, duck and geese, can be used to the hatching of different breeds of eggs by setting different hatching and incubation time. By using intelligent control mode, it can be automatic constant temperature, automatic egg turning and automatic humidifying, with the characteristics of convenient operation, safety and reliability; it can be widely applied to the hatching of egg in various laboratories.

**Technical Specification :**

Hand Drill



**Laboratory Instrument India**