



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Welding Gas Cylinder Trolley

**Product Code :**  
LBNY-0004-930002



**Description :**

Welding Gas Cylinder Trolley

**Technical Specification :**

Trolley for two 10L cylinders (oxygen and acetylene)  
Load Capacity: Two Oxygen Cylinder  
Material: Stainless Steel  
Usage/Application: Industrial  
Color: Blue



**Laboratory Instrument India**