



Email : sales@laboratoryinstrumentindia.com

Product Name :
Digital Multimeter (AC 500V, DC 500V, 200M?)

Product Code :
LBNY-0002-620017



Description :

Digital Multimeter (AC 500V, DC 500V, 200M?)

Technical Specification :

This device gives you accurate reading, helping you measure the electric supply. Well with this product you can measure voltage, current, and resistance of a circuit. The other name for this device is volt-ohm-milliammeter.

This multimeter can be used to measure batteries, household wirings, electric motors and power supplies.

Digital Multimeter, Digital Multimeter Equipments, Digital Multimeter Tools, Digital Multimeter Tool Kits, Digital Multimeter Manufacturers, Digital Multimeter Suppliers from India, China, Kenya



Laboratory Instrument India