



Email : sales@laboratoryinstrumentindia.com

Product Name :
Welding Oxygen Cylinder

Product Code :
LBNY-0004-930003



Description :

Welding Oxygen Cylinder

Technical Specification :

Oxygen cylinder
Usage: Medical Use
Working Pressure: 150 kgf/cm²
Type of Cylinder: Filled
Water capacity: 10L
Working Pressure: 150 bar
Testing Pressure: 250 bar
Cylinder Valve: CGA540
Material: Steel
Protective Cap: Yes



Laboratory Instrument India

