



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Welding Oxygen Cylinder

**Product Code :**  
LBNY-0004-930003



**Description :**

Welding Oxygen Cylinder

**Technical Specification :**

Oxygen cylinder  
Usage: Medical Use  
Working Pressure: 150 kgf/cm<sup>2</sup>  
Type of Cylinder: Filled  
Water capacity: 10L  
Working Pressure: 150 bar  
Testing Pressure: 250 bar  
Cylinder Valve: CGA540  
Material: Steel  
Protective Cap: Yes



**Laboratory Instrument India**

