



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Welding Acetylene Cylinder

**Product Code :**  
LBNY-0004-930004



**Description :**

Welding Acetylene Cylinder

**Technical Specification :**

Acetylene Cylinder  
Water Capacity: 10L  
2.Working Pressure: 150 bar  
Testing Pressure: 250bar  
Cylinder valve: CGA510  
Material: Steel  
Protective Cap: Yes



**Laboratory Instrument India**