



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bio Suit

Product Code :
LBNY-00012-10100052



Description :

The product consists of a 60MM plastic funnel, a support base, a rotatable stainless steel pole, a quantitative suction ball, a mortar, an asbestos net, a goggles, a 6 hole 6 column tube frame, a pair of clamp, 1M silicone tube, a wooden handle tube clip a hoop, belt clip, a plastic spoon, and a cross clamp.

Technical Specification :

Bio Suit



Laboratory Instrument India