



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Just Bench

**Product Code :**  
LBNY-0001-940002



**Description :**

Just Bench

**Technical Specification :**

Just bench  
Simple and economic solutions to the need of production assembly, maintenance and despatch personnel.  
Robust uprights and light grey beam finish.  
Size: 920x1800x750mm  
18mm FSC compliant chipboard working surface.  
Shelves on workstations carry 150kg  
Supplied with metal feet.



**Laboratory Instrument India**