



Email : sales@laboratoryinstrumentindia.com

Product Name :
Locker

Product Code :
LBNY-0001-940005



Description :

Locker

Technical Specification :

Locker

Durable Construction: Made of 16 gauge steel. Includes recessed lift handles with padlock hasp.

Metal locker

Double Tier

Material: Metal

Product Dimensions: 12"D x 36"W x 78"H

Mounting Type Wall Mount

Door Style: Flat Panel

Including Shelves

Shape Rectangular

Base Type Legs

Back Material Type Alloy Steel

Product Dimensions: 12 x 78 x 36 inches

Number of Compartments: 6

Laboratory Instrument India

