



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Locker

**Product Code :**  
LBNY-0001-940005



**Description :**

Locker

**Technical Specification :**

Locker

Durable Construction: Made of 16 gauge steel. Includes recessed lift handles with padlock hasp.

Metal locker

Double Tier

Material: Metal

Product Dimensions: 12"D x 36"W x 78"H

Mounting Type Wall Mount

Door Style: Flat Panel

Including Shelves

Shape Rectangular

Base Type Legs

Back Material Type Alloy Steel

Product Dimensions: 12 x 78 x 36 inches

Number of Compartments: 6

---

# Laboratory Instrument India

