



Email : sales@laboratoryinstrumentindia.com

Product Name :
Handheld LCR Meter

Product Code :
LBNY-0002-620018



Description :

Handheld LCR Meter

Technical Specification :

26-segment analog bar graph; Series/parallel mode (selectable).

Open/short circuit calibration compensation & Smart selection mode for fast selection of components.

5 frequency tests selectable (100Hz/120Hz/1kHz/10kHz/100kHz).

Mini-USB transmission to PC.

Loss factor (D), quality factor (Q), phase angle (?), equivalent series resistance (ESR), equivalent parallel resistance (EPR) and DCR DC resistance measurements.

Inductance (H): 20uH?2000H

Capacitance (F): 200pF?20mF

Resistance (?): 20??200M?

Handheld LCR Meter, Handheld LCR Meter Equipments, Handheld LCR Meter Tools, Handheld LCR Meter Tool Kits, Handheld LCR Meter Manufacturers, Handheld LCR Meter Suppliers from India, China, Kenya



Laboratory Instrument India

