



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Color Polarization Experiment Device For li

**Product Code :**  
LBNY-00014-110007



**Description :**

The circumferential angles of the disc surface and the bottom of the barrel have an angle scale of 0-360 degrees. The instrument consists of two circular polarizers. A mosaic in the bottom of a barrel, another piece is inlaid in the annular disk, such as a disk lid buckle on the drums and can rotate 360 degrees. Outline size 31cm \* 25cm \* 38cm.

**Technical Specification :**

Color Polarization Experiment Device For li



**Laboratory Instrument India**