



Email : sales@laboratoryinstrumentindia.com

Product Name :
Digital Microscope

Product Code :
LBNY-00015-110001



Description :

The item shall have the following specifications.

1. Objective Lens: Glass 4x, 10x, and 40x (retractable) mounted on a revolving objective turret.
2. Stage: includes locked-on stage clips that holds the glass slides; comes with a focus stop to prevent; the objective lens from breaking;
3. Total Magnification (Min.): 40x, 100x, 400x;
4. LCD Screen: 3.5 " LCD screen; displays 5 megapixel images that can be saved to internal memory and to an SD card;
5. Includes minimum of 6 position diaphragm;
6. Camera captures 5 megapixel CMOS images or higher;
7. Built-in minimum of 120 MB Memory; With SD card slot
8. With analog output to connect TV or LCD projector.

Technical Specification :

Digital Microscope



Laboratory Instrument India

