



Email : sales@laboratoryinstrumentindia.com

Product Name :
Digital Microscope

Product Code :
LBNY-00015-110002



Description :

Standard Outfit:

45mm achromatic objective: 4X, 10X, 40(S), 100X (S.Oil);

Wide field plane-scope eyepiece: WF10X; Total magnification:40X-1000X;

Nosepiece?backward quadruple nosepiece;

Double layers mechanical stage: 125x125mm,moving range X-Y: 60x40mm;

?2-?30mm iris diaphragm and ?32 filter; ABBE condenser: N.A.1.25; Coaxial coarse and fine focusing adjustable mechanism: 14mm;

Articulated free binocular head, interpupillary distance: 55-75mm; 30°inclined, 360°rotatable; Built-in digital camera system with a USB2.0 cord and micro-image process software;

Built-in adjustable brightness halogen lamp 6V/20W (or incandescent, or fluorescent, or LED);

Net weight/set: 5kg.

Technical Specification :

Digital Microscope



Laboratory Instrument India

