



Email : sales@laboratoryinstrumentindia.com

Product Name :
Turbidity Meter

Product Code :
LBNY-00016-110002



Description :

Perfect for the plant, field or lab, this simple yet powerful turbidity meter offers an impressive array of features. Use the Infrared (IR) LED light source to comply and reliably measure turbidity with auto-ranging values from 0 to 1000 NTU. Generate an accurate turbidity reading quickly with the Turbidity Meter. This IP67 waterproof meter offers a long battery life averaging about 3000 tests.

Features:

- Auto ranging from 0 to 1000 NTUs
- Uses Nephelometric principle with 90° detection
- Easy to read results from the large, bright LCD display in a variety of work environments
- No mixing or diluting with EPA-approved polymer bead primary standards for faster calibration.

Technical Specification :

Turbidity Meter



Laboratory Instrument India

