



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Turbidimeter

**Product Code :**  
LBNY-00016-110003



**Description :**

Utilize advanced features with the Turbidity Meter, including dual light sources and drinking water/wastewater methods approved for U.S. EPA regulatory reporting. This IP67 waterproof meter records 100 data points and offers an average battery life of 2500 hours intermittent use. Log data into internal memory, or order AQ4CBL cable to send data to printer or download to a computer for easy time- and date-stamped results storage. Keep pace with changing regulatory protocols and use White Light LED to comply with EPA 180.1 requirements or Infrared (IR) LED to comply with ISO 7027 requirements. Select from measure modes including EPA 180.1, ISO Nephelometric, ISO Absorbance, % transmittance, EBC and ASBC. Nephelometric and ratio metric measurements with auto ranging IP67 waterproof portable turbidity meter Dual source to comply with EPA 180.1 and ISO 7027 Readings in the range of 0 to 1000 NTU with a choice of units: NT, FTU, FNU, ASBC or EBC.

**Technical Specification :**

Turbidimeter

---

# Laboratory Instrument India

