



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Turbidity Meter Accessories

**Product Code :**  
LBNY-00016-110004



**Description :**

Ensure your Turbidity Meter has the accessories it needs with our selection of Turbidity Meter Accessories that include turbidity vials, primary turbidity standards, Formazin turbidity standards, and more. Our turbidity meter has advanced features, including dual light sources and drinking water / wastewater methods approved. Accessories for the turbidity meter include primary standards, Formazin turbidity standard, silicon oil and cloth kit, and replacement turbidity vials. Keep pace with changing regulatory protocols; use the White Light LED to comply with EPA 180.1 requirements or Infrared (IR) LED to comply. Select from measure modes including EPA 180.1, Nephelometric, Absorbance, % transmittance, EBC and ASBC. This IP67 waterproof meter records 100 data points and offers an average battery life of 2500 hours intermittent use.

**Technical Specification :**

Turbidity Meter Accessories



**Laboratory Instrument India**

