



Email : sales@laboratoryinstrumentindia.com

Product Name :
Infrared Lab Turbidimeter

Product Code :
LBNY-00016-110006



Description :

Data-logging feature enables user to store up to 1,000 data sets with the date, time, and registration number. Auto-ranging function automatically selects the correct range for turbidity measurement from 0.01 to 1100 NTU. Instrument is power efficient; operates on seven NiCd rechargeable batteries (sold separately) or wall plug adapter.

Technical Specification :

Infrared Lab Turbidimeter



Laboratory Instrument India