



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Turbidimeter Kit

**Product Code :**  
LBNY-00016-110007



**Description :**

The of turbidimeters offers simplistic operation combined with intelligent instrument engineering to provide an unparalleled level of accuracy in turbidity measurement. Featuring the patent-pending sensor technology, the optical system is engineered with dual detectors to deliver a ratio reading which mitigates common measurement stability issues. Reads up to 4000 NTU without the need for dilutions  
Touchscreen interface for ease of operation  
Supplied with AA  
Superior data management - easily save and recall measurement data sets  
Rechargeable battery pack accessory available  
This user interface helps eliminates common frustrations and prevents errors. The data logging capabilities allow you to record the testing location, operator's identification, time and date along with the measurement. The intuitive, touchscreen interface makes it easy to perform procedures and interpret results. With function applications, you can navigate between settings and performing that function efficiently. Stored data can be transferred to a computer via USB.

**Technical Specification :**

Turbidimeter Kit

---

# Laboratory Instrument India

