



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Laboratory Floc / Jar Tester

**Product Code :**  
LBNY-00016-110009



**Description :**

Flocculation Testers (also known as Jar Testers) are designed for a range of applications - such as testing the efficiency of flocculation or precipitation agents  
Digital display of stirring RPM.  
Timer function.  
Date and Time.  
Illuminated back panel.  
Height adjustment option for the stirring blades during operation.  
Continuously variable stirring speed (lab versions only.)

**Technical Specification :**

Laboratory Floc / Jar Tester



**Laboratory Instrument India**