



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Portable Floc Tester

**Product Code :**  
LBNY-00016-1100010



**Description :**

Flocculation Testers (also known as Jar Testers) are designed for a range of applications - such as testing the efficiency of flocculation or precipitation agents  
Illuminated back panel  
Height adjustment option for the stirring blades during operation  
Continuously variable stirring speed (lab versions only).  
Digital display of stirring RPM  
Timer function  
Date and Time.

**Technical Specification :**

Portable Floc Tester



**Laboratory Instrument India**