



Email : sales@laboratoryinstrumentindia.com

Product Name :
Portable Floc Tester

Product Code :
LBNY-00016-1100010



Description :

Flocculation Testers (also known as Jar Testers) are designed for a range of applications - such as testing the efficiency of flocculation or precipitation agents

Illuminated back panel

Height adjustment option for the stirring blades during operation

Continuously variable stirring speed (lab versions only).

Digital display of stirring RPM

Timer function

Date and Time.

Technical Specification :

Portable Floc Tester



Laboratory Instrument India