



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Micro TPW Portable White Light Turbidimeter

**Product Code :**  
LBNY-00016-1100012



**Description :**

Autorange from 0 to 1100 NTU with accuracy of  $\pm 2\%$  of reading or  $\pm 0.01$  NTU  
Water-resistant, ABS-injection molded housing maintains laboratory performance in field operations  
Resolution varies with turbidity range: 0.01 NTU for 0–99 NTU range, 0.1NTU for 100–999.9 NTU range, and 1 NTU for 1000–1100 NTU range  
Light source: tungsten  
Maximum sample size: 10mL  
Other field-friendly features include push-button calibration, four-digit (seven-segment) LCD display, and long life—performs over 1000 tests on a single set of four “AAA” alkaline batteries  
Response time: 14 seconds

**Technical Specification :**

Micro TPW Portable White Light Turbidimeter



**Laboratory Instrument India**