



Email : sales@laboratoryinstrumentindia.com

Product Name :
Accessories for Turbidity Meter

Product Code :
LBNY-00016-1100013



Description :

Generate an accurate turbidity reading quickly with the 0 Turbidity Meter and its accessories. Perfect for the plant, field or lab, this simple yet powerful turbidity meter offers an impressive array of features. This IP67 waterproof meter offers a long battery life averaging about 3000 tests. Use the Infrared (IR) LED light source to comply and reliably measure turbidity with auto-ranging values from 0 to 1000 NTU.

Technical Specification :

Accessories for Turbidity Meter



Laboratory Instrument India