



Email : sales@laboratoryinstrumentindia.com

Product Name :
Turbidity Standard, 200 NTU

Product Code :
LBNY-00016-1100018



Description :

Save time and eliminate errors - no need for dilution

Long shelf life

EPA approved primary turbidity standard

As a primary standard, Turbidity Standards can be used to calibrate or verify any turbidimeter.

All concentrations of Turbidity Standards are verified to meet tight quality control specifications, for reporting purposes and guaranteed to save you time.

Technical Specification :

Turbidity Standard, 200 NTU



Laboratory Instrument India