



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
IR Turbidimeter Kit

**Product Code :**  
LBNY-00016-1100020



**Description :**

Reads up to 4000 NTU without the need for dilutions

Touchscreen interface for ease of operation

The Series of turbidimeters offers simplistic operation combined with intelligent instrument engineering to provide an unparalleled level of accuracy in turbidity measurement.

Superior data management - easily save and recall measurement data sets

Rechargeable battery pack accessory available

Supplied with AA

Featuring the patent-pending Multipath 90° sensor technology, the optical system is engineered with dual detectors to deliver a ratio reading which mitigates common measurement stability issues.

This user interface helps eliminates common frustrations and prevents errors. The data logging capabilities allow you to record the testing location, operator's identification, time and date along with the measurement. Stored data can be transferred to a computer via USB. The intuitive, touchscreen interface makes it easy to perform procedures and interpret results. With function applications, you can navigate between settings and performing that function efficiently.

**Technical Specification :**

IR Turbidimeter Kit

---

# Laboratory Instrument India

