



Email : sales@laboratoryinstrumentindia.com

Product Name :
IV Chemistry Reagents

Product Code :
LBNY-00016-1100022



Description :

The auto-test ampoules contain a pre-measured amount of liquid reagent and easily fill with the correct volume of sample for simple and fast sample and reagent mixing. Each ampoule contains an Auto-ID barcode that is read by the colorimeter and used to automatically identify the test parameter and initiates the test. Minimize the risks of improper measurement and contamination with IV Chemistry Reagents, which are contained in self-filling vials. Use these reagents with the colorimeter to simplify sample preparation and tests.

No. of test: 30

Consistency and reproducibility improved with Aquafast ampoules and an meter
Eliminates measurement errors.

Technical Specification :

IV Chemistry Reagents



Laboratory Instrument India