



Email : sales@laboratoryinstrumentindia.com

Product Name :
PC Application Trainer System

Product Code :
LBNY-0001-950005



Description :

PC Application Trainer System

Technical Specification :

PC Application Trainer System.

Features:

- Easy identification of different parts and components of the system at a glance
- Learn installation, troubleshooting and fault finding
- Study different I/O ports of mother board of PC/AT computer on a PCB
- Easy measurement of voltages and observation of waveforms on test point
- More than 80 test point and 60 switched faults
- Soldering free fault creation and troubleshooting.

Scope of Learning:

- To understand multimedia hardware system
- To observe different cables and connectors
- Learn BIOS configuration
- Hard disk installation and configuration
- To create and delete partitions of hard disk
- General backup techniques
- Diagnosing and fixing hard disk problems

Installing a CD-ROM, CD-R, CD-RW, or DVD Drive.
Configuration of internet
Installing a sound card
Handling I/O address conflicts
Using add new hardware wizard
Troubleshooting of printer
Solving mouse and keyboard problems
Power safety: Line noise protection and battery backup
Installing memory module like RAM and HDD
Learn to disassemble and assemble the computer.



Laboratory Instrument India