



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
IV Chemistry Reagents

**Product Code :**  
LBNY-00016-1100024



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The auto-test ampoules contain a pre-measured amount of liquid reagent and easily fill with the correct volume of sample for simple and fast sample and reagent mixing. Minimize the risks of improper measurement and contamination with IV Chemistry Reagents, which are contained in self-filling vials. Use these reagents with the colorimeter to simplify sample preparation and tests. Each ampoule contains an Auto-ID barcode that is read by the colorimeter and used to automatically identify the test parameter and initiates the test.

Consistency and reproducibility improved with ampoules and an meter.

No. of test: 30

**Technical Specification :**

IV Chemistry Reagents



**Laboratory Instrument India**