



Email : sales@laboratoryinstrumentindia.com

Product Name :
Standard Set B, For IR

Product Code :
LBNY-00016-1100031



Description :

Supplied in ready-to-use sealed vials
Save time and eliminate errors - no need for dilution
All concentrations are verified to meet tight quality control specifications, are EPA compliant for reporting purposes and are guaranteed to save time.
Turbidity Standards are stabilized Formazin standards prepared in ready-to-use concentrations to ensure accurate turbidity measurements.

Technical Specification :

Standard Set B, For IR



Laboratory Instrument India