



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Turbidity Standard, 0.3 NTU

**Product Code :**  
LBNY-00016-1100035



**Description :**

Turbidity Standards are stabilized standards prepared in ready-to-use concentrations to ensure accurate turbidity measurements.

Save time and eliminate errors - no need for dilution

EPA approved primary turbidity standard

Eliminate the opportunity for inaccurate turbidimeter calibrations resulting from the improper dilution or handling of stock solution.

Long shelf life

As a primary standard, Turbidity Standards can be used to calibrate or verify any manufacturer's turbidimeter.

All concentrations of Turbidity Standards are verified to meet tight quality control specifications, are EPA and ISO compliant for reporting purposes and guaranteed to save you time.

**Technical Specification :**

Turbidity Standard, 0.3 NTU



**Laboratory Instrument India**

