



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Replacement Carrying Case

**Product Code :**  
LBNY-00016-1100041



**Description :**

The bottom insert which holds the turbidimeter and key accessories is made of 100% recycled ABS plastic, and can be recycled in the event it becomes damaged and is no longer usable. This set replaces the Carrying Case, upper foam insert and lower tray insert for use with all versions of the turbidimeters. Pieces are not available separately. The case features an easy to grip handle, and the integrated hinges make it robust and field-ready.

**Technical Specification :**

Replacement Carrying Case



**Laboratory Instrument India**