



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Vial Stand for 6 x 24 mm Vials

**Product Code :**  
LBNY-00016-1100043



**Description :**

This Vial Stand offers space for six round cuvettes with a diameter of 24 mm  
While in the stand, the vials are safely stored and cannot slip or tip over - it's the perfect assistant for smooth test procedures!  
Thanks to the transparent acrylic, users always have an unobstructed view of sample development and color changes.

**Technical Specification :**

Vial Stand for 6 x 24 mm Vials



**Laboratory Instrument India**