



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vial Stand for 6 x 24 mm Vials

Product Code :
LBNY-00016-1100043



Description :

This Vial Stand offers space for six round cuvettes with a diameter of 24 mm. While in the stand, the vials are safely stored and cannot slip or tip over - it's the perfect assistant for smooth test procedures! Thanks to the transparent acrylic, users always have an unobstructed view of sample development and color changes.

Technical Specification :

Vial Stand for 6 x 24 mm Vials



Laboratory Instrument India