



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Cable, USB-C

**Product Code :**  
LBNY-00016-1100044



**Description :**

A standard industry USB-C Cable that is able to deliver power to the portable turbidimeters when combined with a power supply.

For Use With (Application): Alternate power supply option for turbidimeters.

For Use With (Equipment): Portable Turbidimeter.

Type: USB-C Cable

Length (Metric): 1 m

**Technical Specification :**

Cable, USB-C



**Laboratory Instrument India**