



Email : sales@laboratoryinstrumentindia.com

Product Name :
Turbidity Meter Kits, Tungsten Version

Product Code :
LBNY-00016-1100046



Description :

Multi-detector optical configuration assures long-term stability and minimizes stray light and color interference.
Easy menu driven operation and large LCD display
Advanced calibration algorithms.
Ideally suited for low-level drinking water applications, mid-range industrial applications, and high-range environmental applications.
500 point data log; stored results can be viewed directly on instrument.

Technical Specification :

Turbidity Meter Kits, Tungsten Version



Laboratory Instrument India