



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Portable Turbidimeters Reusable Glass Cuvettes and Stand

**Product Code :**  
LBNY-00016-1100047



**Description :**

Hard-tempered glass surface resists surface scratching  
Reusable 28 x 70mm glass cuvettes  
Cuvette stand stores up to 11 cuvettes or calibration standards.  
Supplied with threaded light shield caps.

**Technical Specification :**

Portable Turbidimeters Reusable Glass Cuvettes and Stand



**Laboratory Instrument India**