



Email : sales@laboratoryinstrumentindia.com

Product Name :
Portable Turbidimeters Reusable Glass Cuvettes and Stand

Product Code :
LBNY-00016-1100047



Description :

Hard-tempered glass surface resists surface scratching
Reusable 28 x 70mm glass cuvettes
Cuvette stand stores up to 11 cuvettes or calibration standards.
Supplied with threaded light shield caps.

Technical Specification :

Portable Turbidimeters Reusable Glass Cuvettes and Stand



Laboratory Instrument India