



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Densitometer

**Product Code :**  
LBNY-00016-1100051



**Description :**

Densitometer is a compact device that can be used to measure turbidity (in the range of 0.5 to 4.0 McFarland units) of cell suspensions in a variety of life science applications.

This is ideal for determining concentration of cells for fermentation process, detecting the susceptibility of micro-organisms to antibiotics, identifying micro-organisms with various test systems, and measuring optical density at fixed wavelength.

This product features a measurement time of 1 second, a precision of 3%, and small footprint and low-profile design.

**Technical Specification :**

Densitometer



**Laboratory Instrument India**