



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Battery Powered Densitometer, 110V, 75 x 115 x 165mm

**Product Code :**  
LBNY-00016-1100054



**Description :**

The densitometers are designed and factory calibrated to measure turbidity in the range of 0.5 to 4.0 units 0.0 to >7.5 units with a small standard deviation.

A compact benchtop densitometer (suspension turbidity meter) for measuring turbidity of cell suspensions in a variety of life science applications.

If required, it can deliver a wider measurement range (up to 15.00 units), but with a greater standard deviation.

Measurement time: 1 second. Precision:  $\pm 3\%$ .

Measurement range: 0.5 to 4.0 units, 0.0 to >7.5 units.

**Technical Specification :**

Battery Powered Densitometer, 110V, 75 x 115 x 165mm



**Laboratory Instrument India**