



Email : sales@laboratoryinstrumentindia.com

Product Name :
GSC Turbidity Tube

Product Code :
LBNY-00016-1100062



Description :

Fill the tube with water and drain it slowly by limiting flow with the pinch clamp on the hose
Stop draining once the black and white Secchi disk on the bottom of the tube becomes visible
GSC Turbidity Tube measures the cloudiness of a fluid
Take note of the water level using the centimeter graduations on the tube.

Technical Specification :

GSC Turbidity Tube



Laboratory Instrument India