



Email : sales@laboratoryinstrumentindia.com

Product Name :
Replacement Glass Test Tubes for 2020we and 2020wi
Turbidity Meters

Product Code :
LBNY-00016-1100063



www.laboratoryinstrumentindia.com

Description :

For 2020we and 2020wi Turbidity Meters

For Use With (Application): Compact handheld water analysis instrument with microprocessor driven optical system and unique optics configuration

Type: Test Tube

Includes: Screw-on cap

Material: Glass

Technical Specification :

Replacement Glass Test Tubes for 2020we and 2020wi Turbidity Meters



Laboratory Instrument India