



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Portable Turbidimeter

**Product Code :**  
LBNY-00016-1100070



**Description :**

Ideal for regulatory monitoring, process control or field use.  
For fast and easy turbidity measurement in both the field or in the laboratory  
Large easy-to-read LCD screen  
Easy calibration using supplied turbidity standards  
Auto-ranging features selects the proper turbidity range for your sample  
Simple five-key operation  
Sample size: 10-15mL  
Power efficient microcircuit design yields 5000 tests on four "AA" alkaline batteries with an estimated 7-10 year bulb life.

**Technical Specification :**

Portable Turbidimeter



**Laboratory Instrument India**