



Email : sales@laboratoryinstrumentindia.com

Product Name :
Portable Turbidimeter

Product Code :
LBNY-00016-1100070



Description :

Ideal for regulatory monitoring, process control or field use.
For fast and easy turbidity measurement in both the field or in the laboratory
Large easy-to-read LCD screen
Easy calibration using supplied turbidity standards
Auto-ranging features selects the proper turbidity range for your sample
Simple five-key operation
Sample size: 10-15mL
Power efficient microcircuit design yields 5000 tests on four "AA" alkaline batteries with an estimated 7-10 year bulb life.

Technical Specification :

Portable Turbidimeter



Laboratory Instrument India