



Email : sales@laboratoryinstrumentindia.com

Product Name :
Turbidity Std. 1000.0NTU 500ml; Platform: Solution;
Parameter: ; Concentration: 1000 NTU

Product Code :
LBNY-00016-1100074



Description :

Turbidity standards are accepted for calibration purposes. These standards are safe to handle in the lab. No special glassware or difficult lab techniques are required for the tests.

The pre mixed stabilized Formazin turbidity standards are true Formazin dilutions developed for use in any turbidimeter.

Its enhanced stability saves time and labor.; Platform: Solution; Parameter: ; Concentration: 1000 NTU; Precise, Pre-mixed Formazin Solutions; Enhanced Stability Saves Time and Labor; Confidence in results standards have precise concentrations as low as 0.060 NTU for conventional turbidimeters.

Technical Specification :

Turbidity Std. 1000.0NTU 500ml; Platform: Solution; Parameter: ; Concentration: 1000 NTU



Laboratory Instrument India