



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Turbidity Standard, 100FNU/FNRU, 60 mL

**Product Code :**  
LBNY-00016-1100077



**Description :**

Description: Turbidity Standard, 100 FNU/FBRU, 60 mL  
For use with Turbidity Meter, Infrared Version  
Volume (English): 2 oz.  
Volume (Metric): 60 mL  
Turbidity: 100 FNU/FNRU Standard

**Technical Specification :**

Turbidity Standard, 100FNU/FNRU, 60 mL



**Laboratory Instrument India**