



Email : sales@laboratoryinstrumentindia.com

Product Name :
Turbidity Standard, 100FNU/FNRU, 60 mL

Product Code :
LBNY-00016-1100077



Description :

Description: Turbidity Standard, 100 FNU/FBRU, 60 mL
For use with Turbidity Meter, Infrared Version
Volume (English): 2 oz.
Volume (Metric): 60 mL
Turbidity: 100 FNU/FNRU Standard

Technical Specification :

Turbidity Standard, 100FNU/FNRU, 60 mL



Laboratory Instrument India