



Email : sales@laboratoryinstrumentindia.com

Product Name :
Sustainability Solutions Power Form Black 5.6 mil
Sustainable Nitrile Exam Gloves

Product Code :
LBNY-00017-110003



Description :

Protect hands—and the environment—with a sustainable single-use glove that combines performance, durability and tactile sensitivity for critical tasks.
Fully textured for wet and dry grip
Touch-screen compatible
Silicone-free for critical environments
Sustainable Technology
Optimal thickness for durability and tactile sensitivity
Tear-resistant beaded cuff eases donning
1.5 exceeds laboratory standards.

Technical Specification :

Sustainability Solutions Power Form Black 5.6 mil Sustainable Nitrile Exam Gloves



Laboratory Instrument India

